

proBNP antibody [15C4], N-term

Cat. No. GTX39521

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ELISA, Sandwich ELISA
Reactivity	Human

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture : GTX39521, Detection : GTX39517 / GTX39518 / GTX14712 / GTX42656

Not tested in other applications.

Calculated MW	15 kDa. (Note)
Product Note	Human proBNP and NT-proBNP.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	9.20 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic N-terminal proBNP peptide, corresponding to amino acids 61-76 coupled to a carrier protein.
Purification	Protein A purified
Conjugation	Unconjugated

For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).