

proANP antibody [BDI620]

Cat. No. GTX39524

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA
Reactivity	Human

Package
10 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note Synthetic human proANP peptide (a.a. 40-55).

Properties

Form	Liquid
Buffer	50mM Tris
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Synthetic human proANP peptide (a.a. 40-55) poly Lysine conjugated ((C)EPNEEAGAALSPLPEV)
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).