

Pneumocystis carinii (Pneumocystis jiroveci), 65/67kDa antibody [B1236M]

Cat. No. GTX39564

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ICC/IF
Reactivity	Pneumocystis carinii

Package
50 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	1:20

Not tested in other applications.

Product Note Reacts with Pneumocystis carinii cysts and trophozoites. Recognizes 2 antigens at 65 and 67 kDa on Western Blot.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Pneumocystis carinii isolates.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).