

Plasmodium falciparum Histidine-Rich Protein 2 (HRP-2) antibody [B1243M]

Cat. No. GTX39597

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Lateral Flow, Sandwich ELISA
Reactivity	Plasmodium falciparum

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Lateral Flow	Assay dependent
Sandwich ELISA	Assay dependent

Not tested in other applications.

Product Note Reacts with HRP-2 from *P. falciparum*. Does not react with other plasmodium species.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	8.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Native Plasmodium falciparum HRP-2.
Purification	Protein A purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).