

## PKC beta 2 antibody

Cat. No. GTX39622

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, Dot
Reactivity	Rat

Package  
100 µl

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:8,000
Dot	1:1,000

Not tested in other applications.

**Calculated MW** 77 kDa. ( [Note](#) )

**Product Note** In immunoblotting, the product reacts with PKC  $\beta$ 2 (80 kDa) using rat brain extract.

## Properties

Form	Liquid
Buffer	PBS, 1% BSA
Preservative	15mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	synthetic peptide corresponding to the C-terminal region of PKC $\beta$ 2 (660-673).
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).