

# MHC Class II antibody [PF6J-6D]

**Cat. No. GTX39696**

|                     |                |
|---------------------|----------------|
| <b>Host</b>         | Mouse          |
| <b>Clonality</b>    | Monoclonal     |
| <b>Isotype</b>      | IgG2b          |
| <b>Applications</b> | WB, FCM, ELISA |
| <b>Reactivity</b>   | Cat            |

**Package**  
250 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

| Suggested dilution | Recommended dilution |
|--------------------|----------------------|
| WB                 | Assay dependent      |
| FCM                | Assay dependent      |
| ELISA              | Assay dependent      |

Not tested in other applications.

**Product Note** This antibody recognizes an epitope within the 29-31kDa beta chain region of cat MHC class II molecules.

## Properties

|                      |                                                                                                                                                                                                                            |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Form</b>          | Liquid                                                                                                                                                                                                                     |
| <b>Buffer</b>        | PBS                                                                                                                                                                                                                        |
| <b>Preservative</b>  | 0.05% Sodium azide                                                                                                                                                                                                         |
| <b>Storage</b>       | Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles. |
| <b>Concentration</b> | 1.0 mg/ml (Please refer to the vial label for the specific concentration.)                                                                                                                                                 |
| <b>Purification</b>  | Protein A purified<br>From tissue culture supernatant                                                                                                                                                                      |
| <b>Conjugation</b>   | Unconjugated                                                                                                                                                                                                               |

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).