

Myelin basic protein antibody [V/h5]

Cat. No. GTX39736

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, RIA
Reactivity	Human, Rabbit, Sheep, Bovine, Guinea pig

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
RIA	Assay dependent

Not tested in other applications.

Product Note Recognises predominantly residues 45-91.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	N-terminal fragment of human MBP (residues 1-118) produced by BNPS-skatole cleavage of pure human MBP, injected coupled to BSA using carbodiimide.
Purification	Protein A purified
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).