

Marco antibody [ED31]

Cat. No. GTX39744

Host	Rat
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-Fr, FCM
Reactivity	Mouse

References (4)

Package

250 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Assay dependent
IHC-Fr	1/200-1/400
FCM	1/25-1/50

Note : Use 10µl of the suggested working dilution to label 10⁶ cells in 100µl.

Not tested in other applications.

Calculated MW 53 kDa. ([Note](#))

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	CHO cells expressing mouse MARCO.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated




For full product information, images and publications, please visit our [website](#).


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note


Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

DATA IMAGES

GTX39744 IHC-Fr Image


IHC-Fr analysis of mouse lymph node tissue using GTX39744 Marco antibody, green in A and Rat anti Mouse CD8 antibody, clone YTS105.18, red in B. C is the merged image with nuclei counterstained blue using DAPI.


GTX39744 IHC-Fr Image

IHC-Fr analysis of mouse lymph node tissue using GTX39744 Marco antibody [ED31].


GTX39744 IHC-Fr Image

IHC-Fr analysis of mouse lymph node tissue using GTX39744 Marco antibody, green in A and Rat anti Mouse CD8 antibody, clone YTS105.18, red in B. C is the merged image with nuclei counterstained blue using DAPI.


GTX39744 IHC-Fr Image

IHC-Fr analysis of mouse lymph node tissue using GTX39744 Marco antibody [ED31].



For full product information, images and publications, please visit our [website](https://www.genetex.com).