

Leishmania Major Surface Protease gp63 antibody [96/26]

Cat. No. GTX39837

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	ICC/IF, ELISA
Reactivity	Leishmania

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1/500-1/1,000
ELISA	1/2,000

Not tested in other applications.

Product Note This antibody recognizes living and fixed parasites, but will not bind to denatured or improperly folded recombinant gp63.

Properties

Form	Liquid
Buffer	Ascites
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant baculovirus expressed L. major gp63
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).