

## Norovirus Capsid antibody [B114M]

Cat. No. GTX39875

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ELISA, Sandwich ELISA
Reactivity	Norovirus

Package  
500 µg

## Applications

## Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Not tested in other applications.

## Product Note

Recognizes the capsid protein of Norovirus (native and recombinant). Does not cross-react with Rotavirus, Adenovirus or Astrovirus.

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified native Norwalk virus, strain 8FIIa.
Purification	Protein A purified
Conjugation	Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).