

# Norovirus Capsid antibody [B113M]

**Cat. No. GTX39879**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	ELISA, Sandwich ELISA
<b>Reactivity</b>	Norovirus

**Package**  
500 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Not tested in other applications.

### Product Note

Recognizes the capsid protein of Norovirus (native and recombinant). Does not cross react with Rotavirus, Adenovirus or Astrovirus.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	3.30 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Purified native Norwalk virus, strain 8FIIa.
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).