

Trident RIPA Lysis Buffer (Low SDS)

Cat. No. GTX400033

Package
250 ml

PRODUCT

Summary

The ingredients of Trident RIPA lysis buffer (low SDS) contains Tris pH 8.0, NaCl, Triton X-100, sodium deoxycholate, and SDS (0.25%).

PROPERTIES

Form

Liquid

Storage

Store at 2-8°C.

Note

For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).

Date 2024 / 05 / 27 Page 1 of 1