

# Trident RIPA Lysis Buffer (Low SDS)

**Cat. No. GTX400033**

**Package**  
250 ml

## PRODUCT

### Summary

The ingredients of Trident RIPA lysis buffer (low SDS) contains Tris pH 8.0, NaCl, Triton X-100, sodium deoxycholate, and SDS (0.25%).

## Properties

### Form

Liquid

### Storage

Store at 2-8°C.

### Note

For *In vitro* laboratory use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.



For full product information, images and publications, please visit our [website](#).