

## Candida albicans antibody (Biotin)

## Cat. No. GTX40005

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, ELISA, Immunoelectrophoresis, IHC
Reactivity	Candida albicans

Package 1 ml

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent
Immunoelectrophoresis	Assay dependent
IHC	Assay dependent

Not tested in other applications.

**Product Note** It might cross-react with other yeasts.

Properties	
Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	4-5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Candida albicans, type A.
Purification	IgG fraction
Conjugation	Biotin
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 09 Page 1 of 2



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 09 Page 2 of 2