

Streptococcus Enterotoxin A antibody (HRP)

Cat. No. GTX40103

Host	Goat
Clonality	Polyclonal
Isotype	lgG
Applications	ICC/IF, ELISA
Reactivity	Streptococcus

Package 100 μl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:20-1:200
ELISA	1:200-1:1,000

Not tested in other applications.

Product NoteType specific carbohydrate of group A Streptococcus. Does not react with other Strep groups.

Properties	
Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.002% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Streptococci, Group A
Purification	Purified by Affinity chromatography
Conjugation	Horseradish peroxidase(HRP)
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 05 Page 1 of 1

€ 886-3-6208988 📻 886-3-6208989 🐷 infoasia@genetex.com