

# Integrin alpha 3 antibody

**Cat. No. GTX40141**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	IHC-Fr, IHC (Formalin)
<b>Reactivity</b>	Human

**Package**  
100 µg

## Applications

### Application Note

FuncS: Use at an assay dependent dilution. IHC-Fr: Use at an assay dependent dilution. Optimal dilutions/concentrations should be determined by the end user.

**Product Note** Specific to alpha 3 integrin. IA3 is useful for studies on inhibition of cell-cell adhesion and cell-substrate adhesion.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	No preservatives
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
<b>Concentration</b>	1.2 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Cell preparation: Transfected human breast epithelial cell line MTSV1-7, hygromycin resistant
<b>Purification</b>	IgG fraction
<b>Conjugation</b>	Unconjugated

**Note** For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).