

L1CAM antibody

Cat. No. GTX40149

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IP
Reactivity	Human

Package
100 µg

Applications

Application Note

IP: Use at an assay dependent dilution. WB: Use at an assay dependent dilution. Predicted molecular weight: 136 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW	140 kDa. (Note)
Product Note	Binds to tumours of neural and glial origin and does not bind to paediatric or adult brain.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.6 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Tissue / cell preparation. 16 week human foetal brain.
Purification	IgG fraction
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).