

Rad51B antibody

Cat. No. GTX40242

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB
Reactivity	Human

References (1)

Package

100 µg

Applications

Application Note

WB: Use at an assay dependent dilution. Predicted molecular weight: 38 kDa. Optimal dilutions/concentrations should be determined by the end user.

Calculated MW	42 kDa. (Note)
Product Note	Does not recognise Rad51C in total HeLa extracts.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length protein (Human) that has been his-tagged and overexpressed in e.coli.
Purification	IgG fraction The antibody was purified on a talon affinity column under denaturing conditions, followed by gel purification on a SDS-Page.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).