

MHC Class I H-2Kd antibody [B3A6-Do7] (Biotin)

Cat. No. GTX40315

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	FCM, Immunoassay
Reactivity	Mouse

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent
Immunoassay	Assay dependent

Not tested in other applications.

Product Note

Recognizes classical class I antigen, detecting a public specificity shared by both loci K and D.

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Tissue / cell preparation (Mouse lymphocytes).
Purification	Protein A purified
Conjugation	Biotin

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).