

Filamentous Phages M13 / fd / F1 antibody [B62-FE2] (Biotin)

Cat. No. GTX40337

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Applications	WB, ELISA
Reactivity	Bacteriophage M13

Package
125 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1:5,000-1:10,000

Not tested in other applications.

Product Note

Binds to an epitope on pVIII (phage coat protein) covering the N-terminal region of g8pAEGDDPAKAAFDLSLQASAT.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.10 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	fd phages from E. coli F+ strain (JM109).
Purification	Protein A purified From tissue culture supernatant
Conjugation	Biotin

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).