

## MDMA (Ecstasy) antibody [B659M]

Cat. No. GTX40391

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2b
Reactivity	Species independent

Package  
200 µg

## Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	5.20 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Methylenedioxymethamphetamine conjugated to BSA.
Purification	Purified by Column Chromatography
Conjugation	Unconjugated
Note	<p>For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.</p> <p>Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.</p>



For full product information, images and publications, please visit our [website](#).