

Mouse Anti-Human IgA (Fc) CH2 domain antibody [15D6]

Cat. No. GTX40481

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA
Reactivity	Human

Package
1 mg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1/500-1/1,000

Not tested in other applications.

Product Note Clone 15D6 recognizes Human IgA Fc (CH2).

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified Human IgA
Purification	Protein G purified
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTx40481 ELISA Image

Sandwich ELISA analysis of human IgA binding using Mouse anti Human IgA (GTx76197) as a capture reagent and biotinylated Mouse anti Human IgA (GTx40481) as a detection reagent with purified human IgA as antigen. Subsequent detection is by HRP conjugated streptavidin and substrate.



For full product information, images and publications, please visit our [website](https://www.genetex.com).