

Lymphocytic Choriomeningitis virus (LCMV) antibody [M 104]

Cat. No. GTX40495

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ICC/IF, IHC-Fr
Reactivity	Lymphocytic Choriomeningitis Virus

Package
5 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ICC/IF	Neat
IHC-Fr	Assay dependent

Not tested in other applications.

Product Note Reacts with LCMV-infected cells (e.g. Vero, HeLa). Reacts with Human and mouse.

Properties

Form	Liquid
Buffer	Tissue culture supernatant
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Immunogen	LCMV isolated from human HeLa cells.
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).