

LH beta antibody [B212M]

Cat. No. GTX40536

Host	Mouse
Clonality	Monoclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 15 kDa. ([Note](#))

Product Note Does not cross-react with FSH, TSH, or HCG.

Properties

Form	Liquid
Buffer	10mM Tris, 150mM NaCl
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	9.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified natural LH, beta.
Purification	Purified by Ion Exchange Chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).