

## S-arrestin antibody [PDS-1]

Cat. No. GTX40621

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, IHC-Fr
Reactivity	Human, Rat, Bovine, Pig

Package

50 µl

## Applications

## Application Note

## Recommended Starting Dilutions:

For IHC: Use at 1:100

Optimal dilutions should be determined experimentally by the researcher.

Calculated MW 45 kDa. ( [Note](#) )

Product Note S-arrestin in retina from porcine, human, bovine &amp; weakly with rat, and shows heteroditic behaviour

## Properties

Form Liquid

Buffer PBS

Preservative 0.1% Sodium azide

Storage Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.

Concentration Batch dependent (Please refer to the vial label for the specific concentration.)

Immunogen Porcine S-antigen

Purification Protein G

Conjugation Unconjugated

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.

For full product information, images and publications, please visit our [website](#).