

LDL antibody [M-34]

Cat. No. GTX40622

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Reactivity	Species independent

Package

25 µg

Applications

Product Note Blocks the interaction of LDL with human skin fibroblasts.

Properties

Form	Liquid
Buffer	PBS, 0.05% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.05 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified Human LDL
Purification	Purified by affinity chromatography From tissue culture supernatant
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).