

Helicobacter pylori antibody (HRP)

Cat. No. GTX40629

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, IHC-P, ELISA
Reactivity	Helicobacter pylori

Package 100 μΙ

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
IHC-P	Assay dependent
ELISA	Assay dependent
Not tested in other applications	

Not tested in other applications.

Whole cell lysate. Has not been absorbed and may react with related microorganisms. **Product Note**

Properties	
Form	Liquid
Buffer	PBS, 1% BSA
Preservative	0.002% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1-2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	ATCC Strain 43504
Purification	Protein A purified
Conjugation	Horseradish peroxidase(HRP)
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our website.

Date 2025 / 11 / 06 Page 1 of 1