

hCG alpha antibody

Cat. No. GTX40712

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	ELISA
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	Assay dependent

Not tested in other applications.

Product Note

Reacts with alpha chains on human Thyroid Stimulating Hormone (TSH), human Chorionic Gonadotropin(hCG), human Luteinizing Hormone (LH), and human Follicle Stimulating Hormone (FSH).

Properties

Form	Liquid
Buffer	15mM Potassium Phosphate, 0.85% NaCl
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	2.11 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Alpha-hCG
Purification	Purified by affinity chromatography
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).