

Rotavirus VP6 antibody [3C10]

Cat. No. GTX40736

Host	Mouse
Clonality	Monoclonal
Isotype	IgG2a
Applications	WB, ELISA, IHC
Reactivity	Rotavirus

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
IHC	Assay dependent

Not tested in other applications.

Product Note

This antibody recognizes p42 inner capsid protein (VP6). It reacts with bovine and human rotavirus, also reacts with monkey rotavirus (SA-11). This antibody recognizes p42 inner capsid protein. It reacts with bovine and human rotavirus, also reacts with monkey rotavirus (SA-11).

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	Batch dependent (Please refer to the vial label for the specific concentration.)
Immunogen	Purified bovine rotavirus, cell culture adapted strain MR.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).