

Apolipoprotein A1 antibody (HRP)

Cat. No. GTX40784

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA, Immunoassay
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	1:200-1:8,000
Immunoassay	Assay dependent

Not tested in other applications.

Calculated MW 31 kDa. ([Note](#))

Product Note Specific binding to Apo A-I.

Properties

Form	Liquid
Buffer	50mM PBS, 100mM NaCl, 1% BSA
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified human Apo AI
Purification	Purified by Affinity chromatography
Conjugation	Horseradish peroxidase(HRP)



For full product information, images and publications, please visit our [website](#).

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).