

# Insulin antibody [3A6] (Biotin)

**Cat. No. GTX40794**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Reactivity</b>	Human, Bovine, Pig

References ( 1 )

Package

250 µg

## Applications

**Product Note** Cross-reacts with human proinsulin, bovine insulin (30% ) and porcine insulin. No cross-reaction with free C-peptide.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	20mM Tris-HCl, 150mM NaCl
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	0.70 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Purified Human Insulin
<b>Purification</b>	Protein A purified
<b>Conjugation</b>	Biotin

## Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).