

Insulin antibody [3A6] (Biotin)

Cat. No. GTX40794

Host	Mouse
Clonality	Monoclonal
Isotype	lgG1
Reactivity	Human, Bovine, Pig

References (1) Package 250 µg

Applications	
Product Note	Cross-reacts with human proinsulin, bovine insulin (30%) and porcine insulin. No cross-reaction with free C-peptide.
Properties	
Form	Liquid
Buffer	20mM Tris-HCl, 150mM NaCl
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	0.70 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified Human Insulin
Purification	Protein A purified
Conjugation	Biotin
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 19 Page 1 of 1