

Digoxin antibody [3301]

Cat. No. GTX40814

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	EIA, RIA
Reactivity	Species independent

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
EIA	Assay dependent
RIA	Assay dependent

Not tested in other applications.

Product Note GTX40814 shows 100% reactivity with digoxin and digitoxin

Properties

Form	Liquid
Buffer	15mM Potassium Phosphate, 150mM NaCl
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	5.78 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Digoxin- BSA
Purification	Purified by DEAE Chromatography
Conjugation	Unconjugated


Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES

GTX40814 Image


For full product information, images and publications, please visit our [website](#).