

Apolipoprotein A2 antibody (HRP)

Cat. No. GTX40825

Host	Goat
Clonality	Polyclonal
Isotype	IgG
Application	ELISA
Reactivity	Human

Package
100 µl

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ELISA	1:200-1:6,000

Not tested in other applications.

Calculated MW 11 kDa. ([Note](#))

Product Note Specifically binds to human Apo AII.

PROPERTIES

Form	Liquid
Buffer	PBS, 100mM NaCl, 1% BSA
Preservative	0.01% Thimerosal
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Human Apo AII
Purification	Purified by Affinity chromatography
Conjugation	Horseradish peroxidase(HRP)

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).