

IL2 antibody [I2-C7]

Cat. No. GTX40835

Host	Rat
Clonality	Monoclonal
Isotype	IgG
Applications	Neutralizing/Inhibition
Reactivity	Mouse

Package
200 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Neutralizing/Inhibition	Assay dependent

Note : block 50–70% of CTLL proliferative response to “5% Con A sup” as a source for mouse IL-2

Not tested in other applications.

Product Note Does not cross react with rat or human IL-2.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mouse IL-2 purified from Con A activated T cell supernatant.
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).