

pan Endothelin antibody [3G10]

Cat. No. GTX40860

Host	Rat
Clonality	Monoclonal
Isotype	IgG2b
Applications	IHC-P, IHC-Fr, IRMA
Reactivity	Human, Guinea pig

References (1)

Package

100 µg

Applications

Application Note

IHC-Fr: Use at an assay dependent dilution. IRMA: Use at an assay dependent dilution. Is unsuitable for IHC-P. Optimal dilutions/concentrations should be determined by the end user.

Product Note

Recognizes C-terminal region of all endothelins (amino acid residues 15-21) when used as a capture antibody in IRMA assays. Does not cross-react with "big" endothelin.

Properties

Form	Liquid
Buffer	PBS
Preservative	10mM Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Recombinant fragment (Human) conjugated to BSA.
Purification	Protein G purified From hybridoma culture supernatant
Conjugation	Unconjugated
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption. Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).