

## Apolipoprotein E antibody

## Cat. No. GTX40874

Host	Rabbit	
Clonality	Polyclonal	
Isotype	lgG	
Applications	WB, Double immunodiffusion, RID	
Reactivity	Mouse, Rat	

Package 100 μl

## Applications

## **Application Note**

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution		Recommended dilution	
WB		1:500-1:1,000	
Double immunodiffusion		Assay dependent	
RID		1:40	
Not tested in other applications.			
Calculated MW	36 kDa. ( <u>Note</u> )		
Product Note	Species Specific: Extremely low recognition of human Apo E. No data available regarding rat activity.		
Properties			
Form	Liquid		
Buffer	Serum		
Preservative	0.05% Sodium azide		
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.		
Immunogen	Purified Mouse Apolipoprotein E from pooled mouse plasma high and very low density lipoproteins.		
Purification	Unpurified		
Conjugation	Unconjugated		
Note	For laboratory resear	rch use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or in.	
		and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to ure or sequence of the product.	



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 10 Page 1 of 1