

Apolipoprotein E antibody

Cat. No. GTX40874

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Application	WB, Double immunodiffusion, RID
Reactivity	Mouse, Rat

Package
100 µl

APPLICATION

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	1:500-1:1,000
Double immunodiffusion	Assay dependent
RID	1:40

Not tested in other applications.

Calculated MW	36 kDa. (Note)
Product Note	Species Specific: Extremely low recognition of human Apo E. No data available regarding rat activity.

PROPERTIES

Form	Liquid
Buffer	Serum
Preservative	0.05% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Immunogen	Purified Mouse Apolipoprotein E from pooled mouse plasma high and very low density lipoproteins.
Purification	Unpurified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).