

Legionella pneumophila antibody

Cat. No. GTX40943

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	ICC/IF, FCM
Reactivity	Legionella pneumophila

References (3)

Package

1 ml

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	1:800
FCM	Assay dependent

Note : Acetone-fixed is recommended

Not tested in other applications.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	4-5 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	A whole cell preparation of Legionella pneumophila; ATCC #33152
Purification	IgG fraction
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).