

Pancreatic alpha Amylase antibody [503]

Cat. No. GTX40990

Host	Mouse
Clonality	Monoclonal
Isotype	IgM
Applications	Immunoassay
Reactivity	Human

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Immunoassay	Assay dependent
Not tested in other applications.	

Product Note Specifically recognizes human pancreatic amylase. Not tested for reactivity with salivary amylase.

Properties

Form	Liquid
Buffer	PBS
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length native protein (purified): Human Pancreatic Amylase
Purification	Purified IgM
Conjugation	Unconjugated

Note For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).