

# Influenza A virus H1N1 HA (Hemagglutinin) antibody [IVC102]

**Cat. No. GTX40999**

<b>Host</b>	Mouse
<b>Clonality</b>	Monoclonal
<b>Isotype</b>	IgG1
<b>Applications</b>	ICC/IF, FCM, ELISA, Hemagglutinin inhibition
<b>Reactivity</b>	Influenza A virus (H1N1)

References ( 2 )

Package

500 µg

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
FCM	Assay dependent
ELISA	Assay dependent
Hemagglutinin inhibition	Assay dependent

Not tested in other applications.

**Product Note** Influenza virus type A hemagglutinin. Reacts with hemagglutinin of H1 serotype.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	PBS
<b>Preservative</b>	0.09% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	2.60 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Purified influenza virus type A strain H1N1.
<b>Purification</b>	Protein G purified
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).