

MHC Class II antibody [C3A3-Do3] (Biotin)

Cat. No. GTX41007

Host	Mouse
Clonality	Monoclonal
Isotype	IgG3
Applications	FCM
Reactivity	Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
FCM	Assay dependent

Not tested in other applications.

Product Note

GTX41007 detects classical class II antigen on cells from mouse strains carrying Ea ^{d,j,k,p,r,u,v}. The antibody cross-reacts with lymphocytes from rat strains in monomorphic pattern.

Properties

Form	Liquid
Buffer	PBS, 0.1% BSA
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C. DO NOT FREEZE.
Concentration	0.50 mg/ml (Please refer to the vial label for the specific concentration.)
Conjugation	Biotin

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).