

E. coli antibody (HRP)

Cat. No. GTX41110

Host	Goat	
Clonality	Polyclonal	
Isotype	IgG	
Applications	ICC/IF, ELISA	
Reactivity	E. coli	

Package 100 μl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent
Not tested in other applications.	

Product Note

Many "O" and "K" antigenic serotypes of Escherichia coli. Antiserum is not absorbed and does cross-reactwith related Enterobacteriaceae.

Properties	
Form	Liquid
Buffer	PBS, 1% BSA
Preservative	No preservatives
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1-2 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Mixture of E. coli serotypes.
Purification	IgG fraction
Conjugation	Horseradish peroxidase(HRP)
Note	For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.
	Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our <u>website</u>.

Date 2025 / 11 / 04 Page 1 of 1