

Factor VII antibody [RFF-VII/2]

Cat. No. GTX41196

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA, RIA
Reactivity	Human

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
RIA	Assay dependent

Not tested in other applications.

Calculated MW 52 kDa. ([Note](#))

Product Note

This antibody recognizes Factor VII. Clone RFF-VII/2 may be used as a capture antibody in immunoassays for Factor VII in combination with clone RFF-VII/1(GTX41198) as a detection reagent.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.09% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1.0 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified Factor VII.
Purification	Protein G purified From tissue culture supernatant
Conjugation	Unconjugated



For full product information, images and publications, please visit our [website](#).

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Note

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).