

Fetuin A antibody

Cat. No. GTX41201

Host	Sheep
Clonality	Polyclonal
Isotype	IgG
Applications	Double immunodiffusion, Immunoelectrophoresis, RID
Reactivity	Human

Package
100 µl

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Double immunodiffusion	Assay dependent
Immunoelectrophoresis	Assay dependent
RID	Assay dependent

Not tested in other applications.

Properties

Form	Liquid
Buffer	Glycine-buffered saline, 0.1% E-Amino-n-Caproic acid, 0.01% Benzamidine, 1mM EDTA
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	19 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Full length native protein (purified): Human Fetuin A, purified from human serum.
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).