

# Fetuin A antibody

**Cat. No. GTX41201**

<b>Host</b>	Sheep
<b>Clonality</b>	Polyclonal
<b>Isotype</b>	IgG
<b>Applications</b>	Double immunodiffusion, Immunoelectrophoresis, RID
<b>Reactivity</b>	Human

**Package**  
100 µl

## Applications

### Application Note

\*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
Double immunodiffusion	Assay dependent
Immunoelectrophoresis	Assay dependent
RID	Assay dependent

Not tested in other applications.

## Properties

<b>Form</b>	Liquid
<b>Buffer</b>	Glycine-buffered saline, 0.1% E-Amino-n-Caproic acid, 0.01% Benzamidine, 1mM EDTA
<b>Preservative</b>	0.1% Sodium azide
<b>Storage</b>	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
<b>Concentration</b>	19 mg/ml (Please refer to the vial label for the specific concentration.)
<b>Immunogen</b>	Full length native protein (purified): Human Fetuin A, purified from human serum.
<b>Conjugation</b>	Unconjugated

### Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).