

Influenza A virus H1 HA (Hemagglutinin) antibody [B213M]

Cat. No. GTX41227

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	WB, ELISA, Immunoassay
Reactivity	Influenza A virus (H1N1)

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	Assay dependent
ELISA	Assay dependent
Immunoassay	Assay dependent

Not tested in other applications.

Product Note Influenza A Virus Hemagglutinin 1 (H1). Reacts with Hemagglutinin of H1 serotype.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	3.80 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Purified Influenza A/NewCaledonia/20/99 (strain H1N1).
Purification	Protein A purified
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTX41227 WB Image

WB analysis of H1N1 virus (Influenza A/New Caledonia/20/99) and H3N2 virus (Influenza A/Shangdong/9/93) using GTX41227 Influenza A virus H1 HA (Hemagglutinin) antibody [B213M] at 5 μ g/mL.



For full product information, images and publications, please visit our [website](https://www.genetex.com).