

Influenza A virus Nucleoprotein antibody [BDI365]

Cat. No. GTX41247

Host	Mouse
Clonality	Monoclonal
Isotype	IgG1
Applications	ICC/IF, ELISA, Sandwich ELISA
Reactivity	Influenza A virus

Package
500 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
ICC/IF	Assay dependent
ELISA	Assay dependent
Sandwich ELISA	Assay dependent

Note : Capture: GTX41245, Detection: GTX41247/GTX41239.

Not tested in other applications.

Product Note Recognizes Influenza A virus. Specific for nucleoprotein.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.1% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. Store at 4°C.
Concentration	2.70 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	Infected Cell Lysate
Purification	Protein A purified
Conjugation	Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).