

TL1A antibody

Cat. No. GTX41272

Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Applications	WB, ELISA
Reactivity	Human, Mouse, Rat

Package
100 µg

Applications

Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution	Recommended dilution
WB	0.5 - 2 µg/mL
ELISA	Assay dependent

Not tested in other applications.

Calculated MW 28 kDa. ([Note](#))

Product Note This polyclonal antibody will not cross-react with TL1.

Properties

Form	Liquid
Buffer	PBS
Preservative	0.02% Sodium azide
Storage	Store as concentrated solution. Centrifuge briefly prior to opening vial. For short-term storage (1-2 weeks), store at 4°C. For long-term storage, aliquot and store at -20°C or below. Avoid multiple freeze-thaw cycles.
Concentration	1 mg/ml (Please refer to the vial label for the specific concentration.)
Immunogen	TL1A antibody was raised against a peptide corresponding to 14 amino acids near the amino terminus of human TL1A. The immunogen is located within amino acids 60 - 110 of TL1A.
Purification	Purified by DEAE chromatography
Conjugation	Unconjugated

Note


For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).

DATA IMAGES



GTx41272 WB Image

WB analysis of PC-3 cell lysate using GTx41272 TL1A antibody.

Working concentration : (A) 0.5, (B) 1, and (C) 2 μ g/ml



For full product information, images and publications, please visit our [website](https://www.genetex.com).