

Avian Influenza A virus H5N1 HA (Hemagglutinin) blocking peptide

Cat. No. GTX41297-PEP
Applications

Neutralizing /Inhibition

Package

50 µg

PRODUCT
Summary

Avian Influenza A virus HA (Hemagglutinin) peptide is used for blocking the activity of Avian Influenza A virus HA (Hemagglutinin) antibody (catalog number GTX41297).

Applications
Application Note

*Optimal dilutions/concentrations should be determined by the researcher.

Suggested dilution
Recommended dilution

Neutralizing /Inhibition

Assay dependent

Not tested in other applications.

Properties
Form

Liquid

Buffer

PBS, 0.1% BSA

Preservative

0.02% Sodium azide

Storage

Store as concentrated solution. Aliquot and store at -20°C or below. Avoid freeze-thaw cycles.

Concentration

200 µg/ml (Please refer to the vial label for the specific concentration.)

Conjugation

Unconjugated

Note

For laboratory research use only. Not for any clinical, therapeutic, or diagnostic use in humans or animals. Not for animal or human consumption.

Purchasers shall not, and agree not to enable third parties to, analyze, copy, reverse engineer or otherwise attempt to determine the structure or sequence of the product.



For full product information, images and publications, please visit our [website](#).